




Implant Offering Guide

			PRECICE® 2		
	LENGTH	STROKE	8.5 mm	10.7 mm	12.5 mm
ANTEGRADE TIBIA – 10° 	155 mm	50 mm	P8.5-50Q155	P10.7-50Q155	-
	160 mm	30 mm	P8.5-30J160	P10.7-30J160	-
	180 mm	50 mm	P8.5-50J180	P10.7-50J180	-
	195 mm	50 mm	P8.5-50C195	P10.7-50C195	P12.5-50C195
	215 mm	50 mm	P8.5-50C215	P10.7-50C215	P12.5-50C215
	230 mm	50 mm	P8.5-50C230	P10.7-50C230	P12.5-50C230
	245 mm	80 mm	P8.5-80C245	P10.7-80C245	P12.5-80C245
	275 mm	80 mm	P8.5-80C275	P10.7-80C275	P12.5-80C275
	305 mm	80 mm	P8.5-80C305	P10.7-80C305	P12.5-80C305
	335 mm	80 mm	P8.5-80C335	P10.7-80C335	P12.5-80C335
	365 mm	80 mm	-	P10.7-80C365	P12.5-80C365
ANTEGRADE FEMUR TROCHANTERIC – 10° 	165 mm	50 mm	P8.5-50M165	-	-
	170 mm	30 mm	P8.5-30K170	P10.7-30K170	P12.5-30K170
	175 mm	50 mm	-	P10.7-50M175	-
	190 mm	50 mm	P8.5-50K190	P10.7-50K190	P12.5-50K190
	215 mm	50 mm	P8.5-50D215	P10.7-50D215	P12.5-50D215
	230 mm	50 mm	P8.5-50D230	P10.7-50D230	P12.5-50D230
	245 mm	80 mm	P8.5-80D245	P10.7-80D245	P12.5-80D245
	275 mm	80 mm	P8.5-80D275	P10.7-80D275	P12.5-80D275
	305 mm	80 mm	P8.5-80D305	P10.7-80D305	P12.5-80D305
	335 mm	80 mm	P8.5-80D335	P10.7-80D335	P12.5-80D335
	365 mm	80 mm	-	P10.7-80D365	P12.5-80D365
ANTEGRADE FEMUR – PIRIFORMIS STRAIGHT 	215 mm	50 mm	P8.5-50B215	P10.7-50B215	P12.5-50B215
	230 mm	50 mm	P8.5-50B230	P10.7-50B230	P12.5-50B230
	245 mm	80 mm	P8.5-80B245	P10.7-80B245	P12.5-80B245
	275 mm	80 mm	P8.5-80B275	P10.7-80B275	P12.5-80B275
	305 mm	80 mm	P8.5-80B305	P10.7-80B305	P12.5-80B305
	335 mm	80 mm	P8.5-80B335	P10.7-80B335	P12.5-80B335
	365 mm	80 mm	-	P10.7-80B365	P12.5-80B365
RETROGRADE FEMUR – 10° 	165 mm	50 mm	P8.5-50P165	-	-
	170 mm	30 mm	P8.5-30H170	P10.7-30H170	P12.5-30H170
	175 mm	50 mm	-	P10.7-50P175	-
	190 mm	50 mm	P8.5-50H190	P10.7-50H190	P12.5-50H190
	215 mm	50 mm	P8.5-50A215	P10.7-50A215	P12.5-50A215
	230 mm	50 mm	P8.5-50A230	P10.7-50A230	P12.5-50A230
	245 mm	80 mm	P8.5-80A245	P10.7-80A245	P12.5-80A245
	275 mm	80 mm	P8.5-80A275	P10.7-80A275	P12.5-80A275
	305 mm	80 mm	P8.5-80A305	P10.7-80A305	P12.5-80A305
	335 mm	80 mm	P8.5-80A335	P10.7-80A335	P12.5-80A335
	365 mm	80 mm	-	P10.7-80A365	P12.5-80A365
RETROGRADE FEMUR – STRAIGHT 	215 mm	50 mm	P8.5-50E215	P10.7-50E215	P12.5-50E215
	230 mm	50 mm	P8.5-50E230	P10.7-50E230	P12.5-50E230
	245 mm	80 mm	P8.5-80E245	P10.7-80E245	P12.5-80E245
	275 mm	80 mm	P8.5-80E275	P10.7-80E275	P12.5-80E275
	305 mm	80 mm	P8.5-80E305	P10.7-80E305	P12.5-80E305
	335 mm	80 mm	P8.5-80E335	P10.7-80E335	P12.5-80E335
	365 mm	80 mm	-	P10.7-80E365	P12.5-80E365
UNIVERSAL FEMUR – STRAIGHT 	150 mm	50 mm	P8.5-50N150	-	-
	160 mm	50 mm	-	P10.7-50N160	-
	170 mm	30 mm	P8.5-30U170	P10.7-30U170	P12.5-30U170
	190 mm	50 mm	P8.5-50U190	P10.7-50U190	P12.5-50U190

Implant Offering Guide


LOCKING SCREWS

3.5 mm – GREY	LENGTH	PART #
	20 mm	LSB3-020
	25 mm	LSB3-025
	30 mm	LSB3-030
	35 mm	LSB3-035
	40 mm	LSB3-040
	45 mm	LSB3-045
	50 mm	LSB3-050
	55 mm	LSB3-055
	60 mm	LSB3-060

4.0 mm – BLUE	LENGTH	PART #
	20 mm	LSC4-020
	25 mm	LSC4-025
	30 mm	LSC4-030
	35 mm	LSC4-035
	40 mm	LSC4-040
	45 mm	LSC4-045
	50 mm	LSC4-050
	55 mm	LSC4-055
	60 mm	LSC4-060

5.0 mm – GREEN	LENGTH	PART #
	20 mm	LSC5-020
	25 mm	LSC5-025
	30 mm	LSC5-030
	35 mm	LSC5-035
	40 mm	LSC5-040
	45 mm	LSC5-045
	50 mm	LSC5-050
	55 mm	LSC5-055
	60 mm	LSC5-060
	65 mm	LSC5-065
	70 mm	LSC5-070
	75 mm	LSC5-075

END CAP

UNIVERSAL	SIZE	PART #
	Standard	CPA1-000

SIZE GUIDE FOR NAIL DIAMETER

APPLIES TO ALL P2 AND SHORT GEN I NAILS;
MODELS: A, B, C, D, E, J, K, H, AND U (Two Proximal, Two Distal screw holes)

	8.5 mm NAIL		10.7 mm NAIL		12.5 mm NAIL	
	Proximal	Distal	Proximal	Distal	Proximal	Distal
Locking Screw Size	5.0 mm	3.5 mm	5.0 mm	4.0 mm	5.0 mm	5.0 mm

APPLIES TO SHORT GEN II NAILS;
MODELS: Q, M, P, AND N (One Proximal, One Distal screw holes)

	8.5 mm NAIL		10.7 mm NAIL		12.5 mm NAIL	
	Proximal	Distal	Proximal	Distal	Proximal	Distal
Locking Screw Size	5.0 mm	4.0 mm	5.0 mm	4.0 mm	N/A	N/A

For more information about this product, please contact your local sales representative.

101 Enterprise, Suite 100 | Aliso Viejo, CA 92656 Phone: (+1) 949-837-3600 | Fax: (+1) 949-837-3664

CE 0297

©2016, NuVasive, Inc. All rights reserved. NuVasive, and Speed of Innovation are registered trademarks of NuVasive, Inc. NuVasive Specialized Orthopedics, Inc. is a trademark of NuVasive, Inc. PRECICE is a registered trademark of NuVasive Specialized Orthopedics, Inc.

Rx Only.
The PRECICE Intramedullary Limb Lengthening (IMLL) System is composed of an implantable intramedullary nail, locking screws, reusable instruments, and a hand-held External Remote Controller (ERC). The PRECICE nail is a sterile single use device that is surgically implanted using the instruments and locking screws. The ERC is used daily after implantation to non-invasively lengthen or shorten the implant to a prescribed length. The PRECICE System is intended for limb lengthening of the femur and tibia. Contraindications include infection or pathologic conditions of bone such as osteopenia which would impair the ability to securely fix the device, metal allergies and sensitivities, patients whose distance from the surface of the treated limb to the intramedullary canal is greater than 51mm for the 10.7mm diameter implant or greater than 38mm for the 8.5mm diameter implant, patients with an irregular bone diameter that would prevent insertion of the PRECICE nail, patients in which the PRECICE nail would cross joint spaces or open epiphyseal growth plates, patients in which there is an obliterated medullary canal or other conditions that tend to retard healing such as blood supply limitations, peripheral vascular disease or evidence of inadequate vascularity, patients unwilling or incapable of following postoperative care instructions, patients weighing in excess of 114Kg for the 10.7mm diameter implant or weighing in excess of 57Kg for the 8.5mm diameter implant. The implantable device is only to be used by a trained licensed physician. Please refer to the PRECICE IMLL System instructions for use for complete Important Safety Information.