

PediLoc[®]

Locking Proximal
Femur Plate System

*Comprehensive System
for treating pediatric hip
deformity or trauma*

Enhanced fixation

Restoration of mechanical axis


Streamlined, easy to use instruments



Fixation, Alignment, Simplicity


Plate Options and Features

- 3.5mm and 4.5mm plates
- Choice of two offsets for each plate
- 90, 100, 110 degree plates for varus osteotomies
- 130 degree, 6-hole plates for fracture fixation and valgus osteotomies
- 150 degree plates for valgus osteotomies



Screw Options and Features

- Cannulated, converging locking screws for accurate and stable fixation in proximal femur
- 1.6 and 2.0mm guide pins provide control of proximal segment
- Locking and non-locking screws provide options for varying bone quality



OrthoPediatics, Children Are Not Just Small Adults, PediPlates, PediBlade, PediLoc, PediFlex, PLEO, PediNail, Scwire and the Pedi logo are trademarks of OrthoPediatics Corp.

OrthoPediatics, Children Are Not Just Small Adults, PediPlates, PediLoc, Scwire and the Pedi logo are registered trademarks in the United States.

