

PediPlates®

I-Plate and O-Plate

***Indicated for angular
deformity correction***

***Manufactured in
stainless steel***

***Seven different sizes in
2 plate configurations***



Features and Benefits

- Stainless Steel Plates and Screws provide greater strength, resistance to breakage, as well as ease of removal
- Multiple plate sizes and configurations offer greater options intraoperatively
- All screws are self-tapping for easy insertion
- Utilizes 4.5 Stainless Steel Cannulated Screws with multiple screw length options
- Low-profile plates
- The I-Plate more adequately addresses larger and heavier patients by providing additional fixation options
- Dark laser etched line or center hole for guide wire on the plate allows for proper positioning



PediPlates®



The O-Plate is offered in 4 sizes with 2 screw holes per plate.



The I-Plate is offered in 3 sizes with 4 screw holes per plate.



Before correction

During correction

The Adverse Effects and Contra-Indications listed in the IFU are important considerations. These Risks and general surgical Risks should be explained to the patient prior to surgery. †

† IFU, cleaning instructions, and surgical techniques may be obtained by calling OrthoPediatics® Customer Service at 574-268-6379. Read and understand indications, warnings, and adverse effects explained in IFU's prior to use.

OrthoPediatics, Children Are Not Just Small Adults, PediPlates, PediLoc, PediFlex, PLEO, PediNail, Scwire, and the Pedi logo are trademarks of OrthoPediatics Corp.

OrthoPediatics, Children Are Not Just Small Adults, PediPlates, PediLoc, Scwire and the Pedi logo are registered trademarks in the United States. The Pedi logo is a registered Community Trade Mark.



2850 Frontier Drive • Warsaw, IN 46582 • ph: 574.268.6379 or 877.268.6339 • fax: 574.268.6302 • www.OrthoPediatics.com

OrthoPediatics Corp. ©2013 Part# SS-1010-01-01 Rev A