



## PEDIATRIC SPICA TABLE

Often times facilities are limited to using an unstable wood box for casting procedures. OrthoPediatrix has again met the growing, and often over looked, requests of the pediatric orthopedic professional to have a properly engineered device for these procedures.

## FEATURES & SPECIFICATIONS

Often times facilities are limited to using an unstable wood box for casting procedures. OrthoPediatics has again met the growing, and often over looked, requests of the pediatric orthopedic professional to have a properly engineered device for these casting procedures.

- Radiolucent
- Light weight
- Adjustable
- Versatile
- Designed with your patient in mind



COMPONENT	WEIGHT	DIMENSIONS
Patient Platform	9.5 lbs	14" X 14" X 10"
Patient Seat	3.9 lbs	11.75" X 31.75" X 13"
Patient Backboard	.4 lbs	22.5" X 2" X .25"
Arm Boards (2)	.2 lbs	10.75" X 2" X 1"
Assembly	14.2 lbs	36.75" X 26" X 13" MAX